

DISPOSABLE DISTRACTION SCREWS

FOR ACF PROCEDURES

Frontier Devices, a division of Folsom Metal Products manufactures disposable distraction screws with the highest quality Stainless Steel promoting strength and exacting tolerances. Surgeons can have complete confidence in the quality, safety and performance when using them with any vertebral body distractor for anterior cervical fusion.



Self-tapping design for easy three-finger insertion

UNIQUE BENEFITS

- ✓ Compatible with all ACF Systems
- ✓ Unique sharp, self-tapping tip design for three-finger insertion
- ✓ Available in 10 or 20 pack quantities
- ✓ Manufactured of high-grade stainless steel
- ✓ Distraction screws are packaged sterile in single or double packs
- ✓ Available in non-sterile packs

ORDERING INFORMATION

<u>Size</u>	<u>Qty.</u>	<u>Cat. No.</u>	<u>Qty.</u>	<u>Cat. No.</u>
12mm	Ten	301.912	Twenty	301.912S2
14mm	Ten	301.914	Twenty	301.914S2
16mm	Ten	301.916	Twenty	301.916S2

